




Our **Collagen Dental Wound Dressings** are soft resorbable collagen matrices engineered from highly purified Type 1 collagen derived from bovine Achilles tendon. The thickness and pore structure of these devices allow them to absorb fluid and blood at the defect site. In addition, because Collagen Dental Wound Dressings are made of intact collagen fibers, the products have intrinsic hemostatic properties for the control of bleeding.

		
10 mm Øx20 mm	20x40x3 mm	25x75x1 mm
<i>Extraction sites</i> <i>Traumatic wounds</i> <i>Surgical wounds</i>	<i>Periodontal surgical wounds</i> <i>Oral ulcers</i> <i>(non-infected or viral)</i> <i>Denture sores</i>	<i>Extraction sites</i> <i>Traumatic wounds</i> <i>Surgical wounds</i>
Cylindrical shape adapts easily to surgical sites	Thicker sponge shape protects soft tissue	Thinner tape like design for more intricate uses
Can be cut and shaped to size • Individually packed in boxes of 5 units • Single use and supplied sterile		
neo PLUG	neo COTE	neo TAPE



For additional information call your local product specialist or call us toll free at:

1 877 243-6724

www.citagenix.com

