


citagenix

osteodisc[®]

Bone mill for use in grafting procedures

Transform your autologous bone harvest into suitable particle sizes for oral and cranio-maxillofacial grafting procedures.



Simple to operate with ergonomic design

The Osteodisc consists of a non-disposable main unit of medical grade stainless steel and a disposable disc-shaft assembly to ensure the sharpest cutting disc possible.

- **Designed to fit in any sterilizer**
- **Fits surgical air handpieces**
- **Cutting discs available in 500 microns**



“...Its ease of use and exceptional performance make this an indispensable tool for those surgeons requiring particulated autogenous bone ...”

— Jack G. Zosky, D.D.S., F.R.S.D. (C), F.I.C.D. (Toronto, Canada)

For additional information
call your local product specialist
or call us toll free at

1-877-243-6724
www.citagenix.com


citagenix

1105 Autoroute Chomedey
Laval, QC H7W 5J8
www.citagenix.com

