

TITANIUM IMPLANT CARE KITS

Created to maintain the oral health of your implant patients, our solid titanium curettes are specifically designed to remove plaque and calculus without scratching titanium alloy implant components.



For additional information
call your local product specialist
or call us toll free at:

1-877-243-6724

www.citagenix.com

 **Citagenix**

*Innovative and unique solutions
to your surgical needs*



Combo Kit

includes a total of 14 items:

- 3 right titanium curettes in three different angulations**
- 3 left titanium curettes in three different angulations**
- 6 different replacement tips**
- 1 plastic periodontal probe (with 3 tips)**
- 1 Arkansas slip stone**

Extension Kit

includes a total of 8 items :

- 2 right titanium curettes in two different angulations**
- 2 left titanium curettes in two different angulations**
- 4 different replacement tips**

Starter Kit

includes a total of 6 items :

- 1 right titanium curette**
- 1 left titanium curette**
- 2 different replacement tips**
- 1 plastic periodontal probe (with 3 tips)**
- 1 Arkansas slip stone**