


AQUA-PRESS®

Surgical Irrigation System



For additional information
call your local product specialist
or call us toll free at:

1-877-243-6724
www.citagenix.com



*Innovative and unique solutions
to your surgical needs*

AQUA-PRESS® Surgical Irrigation System

The Aqua-Press System provides controlled irrigation of the surgical area through the irrigation syringe.

Aqua-Press system safeguards the sterility of the solution and prevents the risk of air embolism. A standard infusion bag of max. 1000 mL of sterile, physiological (isotonic) saline solution is placed inside the "pressure pouch". The pouch is inflated by the hand pump, which causes the water to flow. Simple and efficient.

Silicone hoses and metal parts can be sterilized up to 135°C

Irrigation during surgical procedures:

Besides improving the visibility of the surgical site; irrigation with a sterile, physiological saline solution may counteract intracellular edema and potential infection.

Sufficient irrigation may accelerate healing, reduce scar tissue formation and minimize post-operative pain.