

HIGH QUALITY. LOW PRICE. BRILLIANT.

Iris™ cordless LED surgical headlamps

Set yourself free with the world's brightest cordless LED surgical headlamp! The IRIS is optimized for the high performance needs of the surgical operating room and utilizes the brightest, whitest LEDs available today.

Each surgical headlamp system includes

Headlamp, rechargeable digitally controlled battery pack and charger, 3 lens clothes, orange light filter, a nylon tote bag and either 2 joysticks (Duo models) or a Solo headlamp adjuster (Solo models).

The digital control on the battery packs allows you to adjust light intensity depending on your needs. These new versatile battery packs can be clipped to your belt or scrubs, put in a pocket or use the web belt provided. For complete freedom the B2X can even be clipped to the headband of your headlamp.

100,000 lux models are recommended for Clinic ENT, Dermatologists, Endodontists, General Surgeons, OB/GYN, Oral, Facial Plastics, Podiatrists, Prosthodontists, Urologists, and Veterinarians or users with 3.5X loupes or higher (higher power loupes require brighter light).

50,000 lux models bring you the best of both worlds.

A completely cordless LED headlamp that provides phenomenally bright, surgical-quality light output, that's comfortable, lightweight and easy to maintain. Recommended for ER physicians, general dentists, dental hygienists, veterinarians or users without loupes, or 3X loupes or lower.

Digitally Controlled Battery Packs

Rechargeable Power Source (Solo models comes with B2X and Duo with B4X)
Lithium-ion, no battery memory with a life of 4 hours per charge

Iris Solo 50: 50,000 lux
Iris Solo 100: 100,000 lux
Iris Duo 100: 100,000 lux
Iris Duo 200: 200,000 lux



Iris is a trademark of and is manufactured by Enova Medical Technologies

PRODUCT MAY DIFFER FROM IMAGE

GP-8043.0 (06/09)

IRIS™

Consistently Bright
LED's last over 50,000 hours with no light reduction ever

•
Easy to Maintain
No fiber optic cable or bulbs to replace

•
Battery Powered
Lithium-ion rechargeable batteries provide 4-8 hours of use

•
Portable
Entire system fits in a drawstring bag

•
Warranty
5 year headlamp warranty

For additional information call your local product specialist or call us toll free at:

1 877 243-6724
www.citagenix.com

Xcitagenix

Innovative and unique solutions to your surgical needs