

PREPARED USING THE PATENTED ALLOWASH® CLEANING PROCESS

RAPTOS

Cancellous and cortical irradiated bone powder

RAPTOS allograft bone offers an osteoconductive scaffold for bone regeneration. Used either alone or as part of a composite bone graft, allograft bone supplies a natural scaffold facilitating the attachment of osteogenic precursor cells.

Raptos allografts are treated using the patented Allowash® technology and are irradiated in their final package for an additional measure of safety.

Allografts	Particle sizes	Available volumes
Cancellous	250 - 1000 µm	0.5, 1, 2 & 5 cc
	125 - 850 µm	0.5, 1, 2 & 5 cc
Cortical	500 - 1000 µm	1, 2 & 5 cc
	1000 - 2000 µm	5 cc
Cortico-cancellous blend	250 - 1000 µm	0.5, 1, 2 & 5 cc



Innovative and unique solutions to your surgical needs



RAPTOS

ALLOGRAFT

Rigorous donor screening, serological testing and formal processing

- *Low dose gamma irradiation*

- *AATB accredited source*

- *Mandatory CTO registered Health Canada*

For additional information you can contact us at:

1 877 243-6724

www.citagenix.com

info@citagenix.com